Actim[®] PROM

BRZI TEST ZA DETEKTOVANJE PREVREMENE RUPTURE FETALNE MEMBRANE

Actim PROM je pouzdan test za brzo i lako detektovanje prevremene ruptуре fetalne membrane (PROM) kod svih pacijenata.

PROM je najčešći identifikovani uzrok prevremenih porođaja i glavni uzrok majčine i fetalne smrtnosti. Ispravna dijagnoza PROM-a je važan vodič tretmana i minimiziranja komplikacija. Ipak, simptomi PROM-a se razlikuju od pacijenta do pacijenta, što čini izazov za dijagnozu.

Sa ACTIM PROM-om, čak i klinički nevidljive ruptуре fetalnih membrana se brzo i tačno detektuju. Za razliku od drugih sličnih testova, Actim Prom daje pouzdane rezultate čak i u slučaju prisustva krvi ili drugih supstanci, čineći test pogodnim za sve žene kod kojih postoji sumnja na PROM. Dobra dijagnoza PROM-a omogućava adekvatan tretman pacijenta, pomaže u izbegavanju nepotrebnog lečenja, i redukuje troškove.

 **Medix
Biochemica**

 www.medixbiochemica.com

KAKO ACTIM PROM RADI

Actim PROM brzi test se zasniva na visoko specifičnim i jedinstvenim monoklonim antitelima koji se vezuju za **insulin-like growth factor binding protein-1 (IGFBP-1)** koji je prisutan u amnionskoj tečnosti tokom trudnoće (slika 2). IGFBP-1 se stvara u ćelijama decidue i glavni je protein u amnionskoj tečnosti. Kada dođe do rupture fetalne membrane amnionska tečnost curi u vaginu, gde se normalno ne nalazi (slika 1). Prisustvo amnionske tečnosti se detektuje sa ACTIM PROM testom.

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Selected references

1. Erdemoglu E and Mungan T. Significance of detecting insulin-like growth factor binding protein-1 in cervicovaginal secretions: Comparison with nitrazine test and amniotic fluid volume assessment. *Acta Obstet Gynecol Scand* (2004) 83: 622–626.
2. Gaucherand P et al. Comparative study of three vaginal markers of the premature rupture of membranes. *Acta Obstet Gynecol Scand* (1997) 76:536–540.
3. Guibourdenche J et al. Rapid detection of insulin-like growth factor-binding protein-1 and foetal fibronectin in cervico-vaginal secretions to diagnose premature membrane rupture. *Ann Clin Biochem* (1999) 36:388–390.
4. Jain K and Morris PG. A clinical study to evaluate the usefulness of the MAST test in diagnosing pre-labour rupture of membranes. *J Obstet Gynaecol* (1998) 18:33–36.
5. Kubota T and Takeuchi H. Evaluation of insulin-like growth factor binding protein-1 as a diagnostic tool for rupture of the membranes. *J Obstet Gynecol Res* (1998) 24: 411–417.
6. Palacio et al. *BMC Pregnancy and Childbirth* 2014, 14:183
7. Ragosch V et al. Insulin like growth factor binding protein 1 (IGFBP-1) und fetales Fibronectin in der Diagnostik eines vorzeitigen Blasensprunges. *GebFra* (1996) 56:1–6.
8. Rutanen E-M et al. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. *Clinica Chimica Acta* (1996) 253: 91–101.
9. Rutanen E-M et al. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROMcheck Membrane immunoassay in the diagnosis of ruptured fetal membranes. *Clinica Chimica Acta* (1993) 214: 73–81.
10. Rutanen E-M. Insulin-like growth factors in obstetrics. *Curr Opin Obstet Gynecol* (2000) 12:163–168.
11. Wathen NC et al. Levels of insulin-like growth factor-binding protein-1 increase rapidly in amniotic fluid from 11 to 16 weeks of pregnancy. *J Endocrinol* (1993) 137:R1–R4.
12. Yang J et al. Vaginal bleeding during pregnancy and preterm birth. *Am J Epidemiol* (2004) 160:118–125.

The full reference list can be found on our website.

KOMBINOVANJE ACTIM PROM SA ACTIM PARTUS

Originalni brzi test za detektovanje prevremene rupture fetalne membrane (PROM), za više samopouzdanja u kreiranju kliničkih odluka

KONTAKT

tel: 011/425 21 21

Besplatan poziv
0800 505051

Ordering information

Actim PROM 20 test kit	30832ETAC
Actim PROM 10 test kit	30831ETAC
Actim PROM 1 test kit	30830ETAC
Actim PROM Controls	30800ETAC
Actim 1ngeni Instrument	19100AC
Actim PROM 1ngeni 10 test kit	30831RETAC